

Drive DeVilbiss® iGo2

Portable Oxygen Concentrator


Preventive Maintenance

Periodic intervals:



Clean and replace the cannula and tubing as instructed by the manufacturer and your equipment provider.

When changing residents:



Clean and disinfect the exterior covers of the device and rechargeable battery between patient use and as needed.



| Specifications: | |
|--------------------------|---|
| Flow Rate: | 1 - 5 pulse settings |
| O ₂ Purity %: | 90% +4/-3% |
| Sound Level: | 37 dBA at setting 2 |
| Battery Life: | 3.5 hours on pulse setting 2 |
| Dimensions: | 3.5"W x 8.6"D x 8.4"H |
| Weight: | 4.95 lbs. |
| Warranty: | 5 years device, 3 years compressor, 1 year sieve beds and battery |

- SmartDose Auto-Adjusting Conserver Technology included
- Settings 1-5 pulse dose
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market